

In door Oxygen Supply System



Application facility



Oxygen concentrator installation cases



Government Office / Public Corporation

- Blue House / NSC / Defense Security Command
- Korea Broadcasting Corporation (KBS 3 locations)
- Asan Water Purification Plant main building
- Administrative agencies, fire department, police station and many others



Apartment / Living space

- Lotte Castle (Yeouido, Seocho, Jamsil)
- GS Construction (Banpo-dong)
- Seoul Forest Prugio Apartment
- Many other apartments




School / Studying facility

- Hankuk University of Foreign Studies Reading Room
- Korea University (Dormitory lounge)
- Gongju Singwan kindergarten
- Education office, many elementary, middle and high schools



General business place

- Lotte Mart, E-Mart, Homeplus nursing rooms
- Marriott Hotel Therapy Room
- Nike Korea
- General office, home 5,000 places, etc.



Sortation	Sauna	Basement	City Office	outside air of the city	Ambient air level	In the woods near the city center	Deep forest
Oxygen level	18~19.3%	20.4% blow	20.5~20.7%	20.8~20.9%	21.0%	21.2%	21.3%
Feeling	Very frustrating	frustrating	Little frustrating	Normal	Refreshing		Very refreshing
Remark	Suffocating	headache	Ventilation required	noxious gas	Optimal oxygen concentration, work efficiency		

Oxygen concentrator models



OXUS 109

- Newest product
- Low noise
- Convenient maintenance

Oxygen flow	10LPM
Oxygen purity	93±3%
Voltage (frequency)	220V(60Hz)
Working temperature	-10~40℃
Power consumption	570W
noise	44dB
size(WxDxH)	390 x 445 x 625mm
weight	40kg



RMC-147C

- Top procurement product in Korea

Oxygen flow	14LPM
Oxygen purity	70%이상
Voltage (frequency)	220V(60Hz)
Working temperature	-10~40℃
Power consumption	570W
noise	60dB
size(WxDxH)	403 x 414 x 622mm
weight	38kg



Controller



Model	CON-RMC1	CON-RMC4
Control method	Contact control	Custom PLC individual power control, 24 hour control
Features	Individual power control and off scheduling	On/Off scheduling
Dimension (WxDxH)	76 x 28 x 121mm	90 x 38 x 120mm

Outlet




Model	DIS-NOZ2	DIS-NOZ7
Features	Display light, humidifier	Ceiling/Wall mount, LED
Dimension (WxDxH)	295 x 100 x 210mm(Acrylic Plate Set)	160 x 38 x 160mm

Indoor air quality controller/outlet




Model	DIS-NOZ8
Features	Oxygen PM, temperature, etc. indication, Auto/manual mode and timer function
Dimension (WxDxH)	660 x 150 x 235mm

Installation layout




History


- 1999. 12 Business Registration – Oxytech
- 2001. 10 Obtaining a medical device manufacturing license
- 2008. 07 IR52 Chang Yong-sil Award (2 times)
- 2013. 03 PPS Selection of Excellent Products for Small and Medium Enterprises
- 2018. 05 New Technology (NET) Certification – Ministry of Industry and Resources
- 2021. 12 Achieving \$2 million in exports in 2021




Certifications and patents_Domestic Patent and Practical Sirian Registration 23 Overseas Patent Registration 9



NET Certi




Korea Certi



Spec Certi



Patent 10-1647017



ISO 13485




ISO9001


Life, Oxygen & Technology




Life, Oxygen and Technology

Address. 1294, Yumyeong-ro, Seorak-myeon, Gapyeong-gun, Gyeonggido, Republic of Korea(ZIP.12469)

 +82-31-477-5252

 sales@oxus.co.kr

 www.oxus.co.kr