

In door Oxygen Supply System



Application facility



Oxygen concentrator installation cases



Government Office / Public Corporation

- Blue House / NSC / Defense Security Command
- Korea Broadcasting Corporation (KBS 3 locations)
- Asan Water Purification Plant main building
- Administrative agencies, fire department, police station and many others



Apartment / Living space

- Lotte Castle (Yeouido, Seocho, Jamsil)
- GS Construction (Banpo-dong)
- Seoul Forest Prugio Apartment
- Many other apartments




School / Studying facility

- Hankuk University of Foreign Studies Reading Room
- Korea University (Dormitory lounge)
- Gongju Singwan kindergarten
- Education office, many elementary, middle and high schools



General business place

- Lotte Mart, E-Mart, Homeplus nursing rooms
- Marriott Hotel Therapy Room
- Nike Korea
- General office, home 5,000 places, etc.



| Sortation | Sauna | Basement | City Office | outside air of the city | Ambient air level | In the woods near the city center | Deep forest |
|--------------|------------------|-------------|----------------------|-------------------------|---|-----------------------------------|-----------------|
| Oxygen level | 18~19.3% | 20.4% blow | 20.5~20.7% | 20.8~20.9% | 21.0% | 21.2% | 21.3% |
| Feeling | Very frustrating | frustrating | Little frustrating | Normal | Refreshing | | Very refreshing |
| Remark | Suffocating | headache | Ventilation required | noxious gas | Optimal oxygen concentration, work efficiency | | |

Oxygen concentrator models



OXUS 109

- Newest product
- Low noise
- Convenient maintenance

| | |
|---------------------|-------------------|
| Oxygen flow | 10LPM |
| Oxygen purity | 93±3% |
| Voltage (frequency) | 220V(50Hz) |
| Working temperature | -10~40℃ |
| Power consumption | 570W |
| noise | 44dB |
| size(WxDxH) | 390 x 445 x 625mm |
| weight | 40kg |



REX-127

- Top procurement product in Korea

| | |
|---------------------|-------------------|
| Oxygen flow | 12LPM |
| Oxygen purity | 70%above |
| Voltage (frequency) | 220V(50Hz) |
| Working temperature | -10~40℃ |
| Power consumption | 520W |
| noise | 60dB |
| size(WxDxH) | 403 x 414 x 622mm |
| weight | 38kg |



Controller



| Model | CON-RMC1 | CON-RMC4 |
|-------------------|---|--|
| Control method | Contact control | Custom PLC individual power control, 24 hour control |
| Features | Individual power control and off scheduling | On/Off scheduling |
| Dimension (WxDxH) | 76 x 28 x 121mm | 90 x 38 x 120mm |

Outlet




| Model | DIS-NOZ2 | DIS-NOZ7 |
|-------------------|--------------------------------------|-------------------------|
| Features | Display light, humidifier | Ceiling/Wall mount, LED |
| Dimension (WxDxH) | 295 x 100 x 210mm(Acrylic Plate Set) | 160 x 38 x 160mm |

Indoor air quality controller/outlet




| Model | DIS-NOZ8 |
|-------------------|--|
| Features | Oxygen PM, temperature, etc. indication, Auto/manual mode and timer function |
| Dimension (WxDxH) | 660 x 150 x 235mm |

Installation layout




History


- 1999. 12 Business Registration – Oxytech
- 2001. 10 Obtaining a medical device manufacturing license
- 2008. 07 IR52 Chang Yong-sil Award (2 times)
- 2013. 03 PPS Selection of Excellent Products for Small and Medium Enterprises
- 2018. 05 New Technology (NET) Certification – Ministry of Industry and Resources
- 2021. 12 Achieving \$2 million in exports in 2021




Certifications and patents_Domestic Patent and Practical Sirian Registration 23 Overseas Patent Registration 9



NET Certi




Korea Certi




Spec Certi



Patent 10-1647017



ISO 13485




ISO9001


Life, Oxygen & Technology




Life, Oxygen and Technology

Address. 1294, Yumyeong-ro, Seorak-myeon, Gapyeong-gun, Gyeonggido, Republic of Korea(ZIP.12469)

 +82-31-477-5252

 sales@oxus.co.kr

 www.oxus.co.kr